RESEARCH

Powered air-purifying respirators used during the SARS-CoV-2 pandemic significantly reduce speech perception

Roxanne Weiss¹, Leon Guchlerner¹, Tobias Weissgerber², Natalie Filmann³, Birgit Haake⁴, Kai Zacharowski⁴, Timo Wolf⁵, Sabine Wicker⁶, Volkhard A. J. Kempf^{7,8}, Sandra Ciesek^{9,10,11}, Timo Stöver¹ and Marc Diensthuber^{1*}

Abstract

Background: Due to the coronavirus disease 2019 (COVID-19) pandemic, interventions in the upper airways are considered high-risk procedures for otolaryngologists and their colleagues. The purpose of this study was to evaluate limitations in hearing and communication when using a powered air-purifying respirator (PAPR) system to protect against severe acute respiratory syndrome coronavirus type 2 (SARS-CoV-2) transmission and to assess the benefit of a headset.

Methods: Acoustic properties of the PAPR system were measured using a head and torso simulator. Audiological tests (tone audiometry, Freiburg speech test, Oldenburg sentence test (OLSA)) were performed in normal-hearing subjects (n = 10) to assess hearing with PAPR. The audiological test setup also included simulation of conditions in which the target speaker used either a PAPR, a filtering face piece (FFP) 3 respirator, or a surgical face mask.

Results: Audiological measurements revealed that sound insulation by the PAPR headtop and noise, generated by the blower-assisted respiratory protection system, resulted in significantly deteriorated hearing thresholds $(4.0 \pm 7.2 \text{ dB hearing level (HL) vs. } 49.2 \pm 11.0 \text{ dB HL}, p < 0.001$) and speech recognition scores in quiet $(100.0 \pm 0.0\% \text{ vs. } 2.5 \pm 4.2\%, p < 0.001$; OLSA: $20.8 \pm 1.8 \text{ dB vs. } 61.0 \pm 3.3 \text{ dB SPL}, p < 0.001$) when compared to hearing without PAPR. Hearing with PAPR was significantly improved when the subjects were equipped with an in-ear headset (p < 0.001). Sound attenuation by FFP3 respirators and surgical face masks had no clinically relevant impact on speech perception.

Conclusions: The PAPR system evaluated here can be considered for high-risk procedures in SARS-CoV-2-positive patients, provided that hearing and communication of the surgical team are optimized by the additional use of a headset.

Keywords: COVID-19, SARS-CoV-2, Personal protective equipment, Powered air-purifying respirator, PAPR, FFP3 respirator

Full list of author information is available at the end of the article

© The Author(s). 2021 **Open Access** This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons licence, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons licence, and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this licence, visit http://creativecommons.org/licenses/by/4.0/. The Creative Commons Public Domain Dedication waiver (http://creativecommons.org/publicdomain/zero/1.0/) applies to the data made available in this article, unless otherwise stated in a credit line to the data.

Weiss et al. Journal of Occupational Medicine and Toxicology (2021) 16:43 https://doi.org/10.1186/s12995-021-00334-y

> Check for updates

Journal of Occupational Medicine and Toxicology



^{*} Correspondence: marc.diensthuber@kgu.de

¹Department of Otorhinolaryngology, Head and Neck Surgery, University Hospital Frankfurt, Goethe University, Theodor-Stern-Kai 7, 60590 Frankfurt/M, Germany

A local outbreak of a previously unknown disease, COVID-19, caused by infection with a novel coronavirus, SARS-CoV-2 [1], in Wuhan, Hubai Province, China, rapidly developed into a global epidemic in early 2020 and is currently posing major challenges to the world's healthcare systems.

This remarkably variable disease shows a broad spectrum of clinical manifestations, ranging from completely asymptomatic patients [2] to rapidly progressive courses with lethal outcome despite intensive care treatment [3]. According to current data, the main route of human-to-human transmission of the pathogen primarily occurs via respiratory droplets from infectious individuals or, less frequently, as a result of direct contact with SARS-CoV-2-contaminated surfaces [4].

The experience of the past year has shown that medical personnel in the operating room are exposed to a considerable risk of infection when treating SARS-CoV-2-positive patients [5]. Since high SARS-CoV-2 viral loads can be detected in the upper respiratory tract of COVID-19 patients [6], otolaryngologists, whose activities are focused on this area, are among the most exposed specialties. A particular high risk of virus transmission appears to emanate from aerosolproducing interventions, such as endonasal skull base surgery using high-speed drill [7] or tracheostomy [8, 9]. Adequate personal protective equipment (PPE) is required not only for operations on COVID-19 patients but also for emergency interventions, e.g. in case of acute respiratory distress or life-threatening bleeding in the upper airways. This is of particular importance when the SARS-CoV-2 test result is unknown or delayed.

For these high-risk interventions N95 respirators (which correspond to the European respirator standards FFP2/3) and goggles have been suggested to protect the surgeons [10]. However, there is some evidence that higher-level PPE in the form of powered air-purifying respirators (PAPRs) may be more effective in safely preventing infection of the surgeon, especially in aerosol-producing, high-risk procedures [11, 12]. As a consequence, the use of PAPRs is now included in some PPE recommendations for interventions during the COVID-19 pandemic [13–15].

In routine clinical practice, we have noticed that hearing and communication are significantly impaired when PAPRs are used. There is a small number of studies providing data which confirm this subjective impression [16, 17]. The purpose of this study was to provide an indepth analysis of acoustic properties of a PAPR which could potentially lead to difficulties in communication. We also tested the impact on hearing and communication when two communication partners used PAPRs to simulate the situation between surgeon and assistant during surgery. The audiological data measured when using PAPRs was compared with those measured when wearing FFP3 respirators or surgical face masks. Moreover, we evaluated the utility of a headset when using a PAPR system as an aid to overcome audiological limitations.

Material and methods

Powered air-purifying respirator (PAPR)

The Tornado T9 respirator headtop (Scott Health and Safety Ltd., Skelmersdale, United Kingdom) (Fig. 1A) with connection to a respiratory blower is a CE-certified device that belongs to the category of PAPRs. The headtop completely encloses the head and is made of translucent polyurethrane, which provides all-round view. It was used in conjunction with a blower-assisted respiratory protection system (PM Proflow 2 SC, PM Atemschutz GmbH, Mönchengladbach, Germany) equipped with a CE-certified PM breathing protection filter CF 32 A2B2E2K2-P3 RD / CF 32 ABEK-P R SL (PM Atemschutz GmbH) suitable to protect against bacteria and viruses, gases and vapours as well as radioactive and highly toxic particles (Fig. 1B). Figure 1C shows the use of the PAPR in combination with an inear headset connected to a cordless telephone.

Measurements of the acoustical properties of the PAPR system

Acoustical properties (sound attenuation and fan noise) of the PAPR system were measured in an anechoic chamber with the system worn on a head and torso simulator (Type 4100, Brüel & Kjær, Nærum, Denmark) connected to a measurement amplifier (Type Nexus, Brüel & Kjær) (Fig. 2A).

The noise level inside the headtop (at the position of the eardrum) generated by the fan of the blower-assisted respiratory protection system was measured in onethird-octave bands. To assess the sound insulation of the headtop, the transfer function from a loudspeaker (C5 tiny, KS Digital, Saarbrücken, Germany) to the left ear of the head and torso simulator was measured with and without the PAPR system. Analysis was conducted with ARTA software (ARTALABS, Kastel Luksic, Croatia).

Participants and audiometric testing

Audiometric testing was performed in normal-hearing subjects (n = 10 health care professionals from our institution; n = 5 males, n = 5 females; age 28.6 ± 4.3 years) who were recruited to participate in the study on a voluntary basis. The participants had no previous regular experience in the use of PAPRs and received detailed instruction in the handling and use of this PPE before the audiometric measurements were performed. All tests were carried out in a sound-isolated booth with regularly





calibrated audiometers and in free-field conditions with the PAPR system (headtop and activated respiratory protection system) and without the PAPR system. To overcome potential hearing deficits while wearing a PAPR, tests were additionally conducted in a simulated communication situation via wireless telephones. In this test condition, the subject used earphones (M2–022 earphone, M2-TEC, China) connected to a wireless telephone device (8242 DECT Handset, Alcatel-Lucent Enterprise, Colombes, France). A second wired microphone (M2–022 earphone, M2-TEC, China) was placed 10 cm (tone audiometry and Freiburg number and monosyllable tests) or 15 cm (speech reception thresholds, German matrix test) in front of the loudspeaker and connected to a second identical wireless telephone.

Warble tone hearing thresholds in free-field with and without PAPR system and with wireless telephone device were measured for the test frequencies 0.125/0.25/0.5/0.75/1/1.5/2/3/4/6/8 kHz. Speech perception scores in quiet were measured using the Freiburg number and monosyllable tests [18] at sound pressure levels of 65 dB SPL corresponding to a medium loud conversational level. Speech was presented from frontal direction at 0°.

Speech reception thresholds (SRTs, sound pressure level of speech with 50% correct word perception) in quiet were measured with the German matrix test (Oldenburg sentence test, OLSA) [19] in an adaptive procedure. Initial sound pressure level of the speech signal was 50 dB SPL. The test was performed in a closed-set mode. The OLSA test was performed without and with PAPR system (headtop and activated respiratory protection system). Additionally, the frequency characteristics of the target speaker were modified to simulate a real conversation, where the person speaking in the operating room would also be wearing PPE. This was realized by convoluting the speech signal of the OLSA test with the frequency response (i.e. dampening function) of three different types of PPE: 3-layer surgical face mask (KF-B P01, Kingfa Science &Tech, China), FFP3 respirator (REF 35100, FarStar Medical, Germany) and PAPR. The convolution was done with the software Equalizer APO (https://sourceforge.net/projects/equalizerapo/). The dampening functions of the 3-layer surgical face mask and the FFP3 respirator are shown in Supplemental Fig. S1. The experimental conditions of the audiometric tests are summarized in Tables 1 and 2.

Statistics

Fig. 2 Measurement of the acoustical properties of the powered air-purifying respirator (PAPR) system. A Illustration of the PAPR system on a head and torso simulator used for measurements. B One-third octave band noise levels (dB SPL) of the PAPR system (black line) and the noise floor of the anechoic chamber (grey line). C Frequency spectrum of the sound insulation (dBr) of the PAPR calculated as difference between the

Quantitative data are given as mean ± standard deviation (and median). Graphical presentation of the data was performed using GraphPad Prism 8 (GraphPad Software, San Diego, USA) and ARTA software (ARTALABS). Statistical data analysis was performed with R (version 4.0.4, R Foundation for Statistical Computing, Vienna, Austria) using linear mixed-effects models. Thereby, a maximum likelihood approach was used for fitting and the AIC (Akaike Information Criterion) was utilized for model selection of fixed and random effects with up to two-fold interactions as proposed by Seedorff et al. [20]. As all data arise from randomized block designs, where each block corresponds to one subject, subject-specific random effects were included to account for this correlation.



Audiological test	Condition
	Listener (subject)
Warble tone audiometry	without PAPR
	with PAPR
	with PAPR and headset
Freiburg number and monosyllable test	without PAPR
	with PAPR
	with PAPR and headset
German matrix test (Oldenburg sentence test, OLSA)	without PAPR
	with PAPR
	with PAPR and headset

transfer functions to the position of the left eardrum with and without PAPR system

(See figure on previous page.)

Audiological test	Condition		
	Listener (subject)	Speaker	
German matrix test (Oldenburg sentence test, OLSA)	SFM	SFM	
	FFP3	FFP3	
	PAPR	PAPR	

Table 2 Audiometric evaluation in subjects using PAPR, FFP3 respirator or surgical face mask with additional frequency response modification of the speaker (manipulation of listener and speaker)

SFM Surgical face mask, FFP3 FFP3 respirator, PAPR Powered air-purifying respirator

Audiometric data were analyzed using the following three models:

- 1. In order to analyze the effect of PAPR (and headset) on hearing thresholds, the conditions "without PAPR", "with PAPR" and "with PAPR and headset", frequency and the interaction between condition and frequency were included as fixed effects. The subject specific effect and the condition were included as random effects.
- 2. The speech recognition scores were analyzed including the condition (without PAPR, with PAPR, with PAPR and headset) as fixed effect and the subject specific effect as random effect.
- 3. To evaluate the influence of PPE on SRT the variable condition (without PAPR, with PAPR, with PAPR and headset, S: SFM / L: SFM, S: FFP3 / L: FFP3, S: PAPR / L: PAPR) was included as fixed and random effect and the subject specific effect as random effect. Hereby we considered the following two-step approach for the pairwise comparisons, as a difference of less than 2 dB is not clinically relevant: first, a test of equivalence for an equivalence margin of ±2 dB was carried out. Then, in case equivalence was not fulfilled, a test of difference was conducted.

For model building, the default "treatment contrasts" of the R language was used, which corresponds to dummy coding. Tukey's all-pair comparisons were used for post-hoc-tests. Post-hoc *p*-values were adjusted with the single-step method. *P*-values ≤ 0.05 were considered statistically significant.

Results

Acoustical properties of the PAPR system

Sound pressure level of the noise generated by the fan of the blower-assisted respiratory protection system is shown in Fig. 2B (black line). The noise floor of the anechoic chamber is also shown (light grey line). One-third octave band noise levels ranged from 32 to 56 dB SPL with highest levels in the frequency region of 250–500 Hz and 1-2 kHz. Sound insulation of the headtop is shown in Fig. 2C. For frequencies up to 400 Hz the headtop can be considered as acoustically transparent. For frequencies higher than 400 Hz an increased sound insulation could be observed with increasing frequency. In the frequency region of 3–5 kHz, which is important for the intelligibility of speech, sound insulation was in the range of 20 dB.

Audiological data

Warble tone audiometry

Hearing thresholds in free-field audiometry without PAPR, with PAPR and with PAPR and headset are shown in Fig. 3A. Average free-field audiometry hearing thresholds for frequencies between 0.125 and 8 kHz without and with the PAPR system were 4.0 ± 7.2 dB HL vs. 49.2 ± 11.0 dB HL (median: 5.0 dB HL vs. 50.0 dB HL). A significant deterioration of the hearing thresholds caused by the sound attenuation of the head top and the fan noise of the activated respiratory protection system was detected (p < 0.001). When the PAPRs were used with a headset a significant improvement of the hearing thresholds was measured for the frequencies between 0.5 and 2 kHz (42.1 ± 5.4 dB HL, median 40.0 dB HL; p < 0.001).

Speech audiometry

Speech recognition scores without PAPR, with PAPR and with PAPR and headset are shown in Fig. 3B. Average number recognition score was $100 \pm 0.0\%$ (median: 100.0%) without PAPR at 65 dB SPL and significantly decreased to $45.0 \pm 15.8\%$ (median: 40.0%) at 65 dB SPL with the PAPR (p < 0.001). Likewise, the average monosyllabic word recognition score without the PAPR was $100.0 \pm 0.0\%$ (median: 100.0%) and was significantly reduced when the PAPR was used ($2.5 \pm 4.2\%$, median: 0.0%, p < 0.001). When using the PAPR system equipped with the in-ear headset, data showed a significant improvement of speech perception to $95.0 \pm 7.1\%$ (median: 100.0%, p < 0.001) for numbers and $52.5 \pm 19.3\%$ (median: 50.0%, p < 0.001) for monosyllabic words.

SRTs obtained with the OLSA for the test conditions "without PAPR", "with PAPR", "with PAPR and headset" and for the simulations of a talker using either surgical face mask, FFP3 respirator or PAPR are shown in Fig.



(See figure on previous page.)

Fig. 3 Impact of the powered air-purifying respirator (PAPR) on hearing performance in normal-hearing subjects (n = 10). **A** Boxplots of hearing thresholds (dB HL) with PAPR (white), without PAPR (grey) and with PAPR + headset (plaid) using tone audiometry. The vertical bar (***) applies for "without PAPR" vs. "with PAPR" and for "without PAPR" vs. "with PAPR and headset". **B** Boxplots of speech recognition scores (%) for numbers and monosyllabic words (Freiburg speech test) without PAPR (black), with PAPR (grey) and with PAPR + headset (white) at 65 dB. **C** Boxplots of speech reception thresholds (SRTs) in dB SPL with PAPR (straight striped), without PAPR (white) and with PAPR and headset (plaid) as well as with attenuation of PAPR (oblique plaid), with attenuation of a surgical face mask (oblique striped) and with attenuation of an FFP3 respirator (grey) obtained with the German matrix test. Box plots show minima, maxima, interquartile range and median. Significant differences not indicated in the graph: without PAPR vs. with PAPR vs. with PAPR and headset, without PAPR vs. S: PAPR / L: PAPR, S: SFM / L: SFM vs. with PAPR, set SFM vs. with PAPR, set SFM vs. with PAPR and headset, S: FFP3 / L: FFP3 vs. with PAPR and headset, S: FFP3 / L: FFP3 vs. with PAPR and headset, S: FFP3 / L: FFP3 vs. S: PAPR / L: PAPR, S: FFP3 / L: FFP3 vs. with PAPR and headset, S: FFP3 / L: FFP3 vs. With PAPR (p < 0.001, respectively). S = speaker, L = listener, SFM = surgical face mask, FFP3 = FFP3 respirator. *p < 0.05; ***p < 0.001

3C. Mean SRTs significantly increased (i.e. deteriorated) from $20.8 \pm 1.8 \text{ dB}$ SPL (median: 21.5 dB SPL) to $61.0 \pm 3.3 \text{ dB}$ SPL (median: 62.0 dB SPL) when the PAPR was used (n = 10, p < 0.001). When the PAPR was equipped with a headset, significantly lower SRTs were observed ($56.2 \pm 4.7 \text{ dB}$ SPL, median: 55.8 dB SPL; p < 0.05).

OLSA results from measurements simulating the effects of different types of PPE used by the speaker showed that SRTs were not affected by the use of either a 3-layer surgical face mask or an FFP3 respirator, as for the comparisons of these conditions ("without PAPR vs. S: SFM / L: SFM"; "without PAPR" vs. "S: FFP3 / L: FFP3"; "S: SFM / L: SFM" vs. "S: FFP3 / L: FFP3") the limits of the simultaneous 95%-confidence intervals for the differences were below 2 dB. Hence, it can be assumed that these conditions are equivalent.

When both speaker and listener (i.e. subject) used a PAPR (condition "S: PAPR / L: PAPR"), mean SRTs significantly deteriorated by 47.7 ± 10.0 dB SPL in comparison to the situation without PAPR (p < 0.001).

Discussion

Since SARS-CoV-2 appears to be transmitted mainly via aerosol particles [21] and droplets [22], interventions in the upper airways, especially if they produce aerosols, are associated with a particularly high risk of SARS-CoV-2 infection for the surgeon as well as the surrounding team and may require enhanced PPE [23, 24]. Consequently, the most exposed specialties are otorhinolaryngology, anesthesiology and dentistry, which together account for 12% of COVID-19-related deaths among physicians [25]. Recent data show that at a university department of otorhinolaryngology in times of the COVID-19 pandemic, approximately 1 in 200 patients were found to be SARS-CoV-2-positive by polymerase chain reaction (PCR) testing [26]. In the field of otolaryngology, in particular the frequently performed routine interventions such as tracheostomy and surgical procedures with powered devices including functional endoscopic sinus surgery and mastoidectomy are considered high-risk procedures [13, 27].

For such interventions, the use of PAPRs has been proposed to ensure enhanced safety for the surgeons involved [13]. The successful use of PAPRs for viral infection control has already been reported in the management of severe acute respiratory syndrome (SARS) patients [28]. Our university hospital also has substantial experience in the use of PAPRs from the care and treatment of highly infectious Ebola patients [29]. Given the increasing incidences of new SARS-CoV-2 variants of concern in Europe [30] and limited availability of vaccines, but also the risk of further new respiratory infectious diseases with pandemic potential in the future, enhanced PPE, such as PAPRs, will continue to play a crucial role in protecting health care professionals from viral transmission.

The three-dimensional protection of the head and neck from splashes and aerosols and the superior filtering capacity of PAPRs compared to N95 respirators and surgical face masks are the defining qualities of PAPRs [11, 24]. Several drawbacks and limitations exist for N95 respirators such as suboptimal respirator fit due to facial hair [31] or significantly increased breathing resistance [32] which must be considered in case of long surgeries in infectious COVID-19 patients. In these cases, PAPRs can be a suitable and valuable alternative.

However, we found in clinical routine practice that hearing and thus communication is significantly compromised when the PAPR is worn and the respiratory protection system has been activated. This limitation is especially critical in case of emergency interventions, where rapid action and optimal communication are key.

To further explore these acoustic difficulties, we conducted hearing tests on normal-hearing subjects with and without PAPR. Tone audiometry revealed that the hearing thresholds were significantly deteriorated when the subjects wore the PAPRs with the blower-assisted respiratory protection system activated. Likewise, speech perception was significantly reduced under these conditions. This specific limitation of PAPRs has been reported previously [11, 16, 17] and has to be considered when this PPE is used in clinical routine. The underlying cause of these results became obvious when we

characterized the acoustic properties of the PAPR. Firstly, the headtop of the PAPR causes substantial sound insulation. Secondly, the activated respiratory protection system produces noise with sound pressure levels up to 56 dB. Both factors together result in the detected restrictions in hearing and communication. Kempfle et al. [16] showed a mean deterioration of hearing thresholds of about 40 dB in a different PAPR system. In our study, the measured hearing thresholds were even 10 dB worse. The obtained mean reduction in word recognition scores of about 50% was comparable with the results reported in the study by Kempfle et al. [16]. Palmiero et al. [17] performed measurements of speech transmission index (STI) for different types of PPE and also reported a significantly lower (i.e. worse) STI for PAPR systems compared with N95 respirators and 3layer surgical face masks.

To compensate for these limitations in hearing, the benefit of in-ear headsets was tested. In tone audiometry, a significant improvement of hearing thresholds of about 10 dB was found in the frequency range between 0.5 and 2 kHz. On the one hand, sound insulation of the PAPR in this frequency range was overcome by the use of a headset. On the other hand, hearing thresholds could not be further improved because of the high fan noise level between 0.25 and 2 kHz. For lower and higher frequencies, no benefit was found by using the headset. This can be explained by the limited transmission of frequencies lower than 500 Hz and higher than 4 kHz via the telephone. Our data (number and word recognition results in quiet) also clearly show that hearing with PAPR is significantly improved when an in-ear headset is used. With PAPR mean word recognition score was almost 0% and improved to about 50% at a conversational level of 65 dB SPL. Taking redundancy of conversational speech into account it could be assumed that a sentence recognition score of 100% is achievable by using a headset. However, listening effort in such a communication situation is still highly increased compared to a conversation using other types of PPE (FFP3 respirator, 3-layer surgical face mask).

This finding was also reflected in the SRT measurements. It was shown that using a PAPR led to a significant deterioration of SRT (40 dB) whereas no clinically relevant impact of FFP3 respirators and 3-layer surgical face masks on SRT was found. In a test condition in which listener and speaker used a PAPR, SRTs increased even more. In the test condition "with PAPR and headset" a significant improvement of 5 dB in SRTs was found. Given that the slope of the OLSA discrimination function is 11.3%/dB, this corresponds to an improvement in speech perception of 56.5%.

In the present study, no clinically relevant impact of surgical face masks or FFP3 respirators on SRT in quiet was found, although the measured sound attenuation in the frequency range between 2 and 4 kHz is up to 10 dB (Supplemental Fig. S1). Other studies also found no deterioration in speech perception for surgical face masks in either quiet [33] or for high signal-to-noise ratios [34]. A subjective decline in communication ability is frequently reported in clinical routine practice when using an FFP3 respirator in comparison to a 3-layer surgical face mask. However, only speech tests in quiet in a sound isolated booth were conducted. It could be hypothesized that differences in speech perception between FFP3 respirator and 3-layer surgical face mask are revealed by conducting a speech test in noise. The major differences in the dampening function between FFP3 respirator and 3-layer surgical face mask are present for frequencies higher than 5 kHz and thus could potentially reduce the recognition of sibilants in noise. Reduced speech perception in noise was shown by Brown et al. [33] for surgical face masks and by Toscano et al. [34] for N95 respirators (dampening function in the range of FFP3 respirator) in test conditions in noise with lower signal-to-noise ratio. In addition, our tests do not measure the potential negative effects of FFP3 respirators on articulation due to the tight fit on the jaw and the lips which could also affect communication ability.

Conclusions

In summary, we present a PAPR system that appears suitable for use in aerosol-generating procedures. However, the use of such a PPE is accompanied by limitations in hearing and communication. If these difficulties are overcome by providing health care professionals with headsets, PAPRs represent a PPE with a very high level of protection, which can be recommended for high-risk interventions during the SARS-CoV-2 pandemic.

Since PAPR systems can be used not only to protect health care professionals from SARS-CoV-2 infection but also from various other highly contagious respiratory diseases, it can be anticipated that our findings will also be relevant for future pandemics.

Abbreviations

COVID-19: Coronavirus disease 2019; FFP: Filtering face piece; HL: Hearing level; OLSA: Oldenburg sentence test; PAPR: Powered air-purifying respirator; PCR: Polymerase chain reaction; PPE: Personal protective equipment; SARS: Severe acute respiratory syndrome; SARS-CoV-2: Severe acute respiratory syndrome coronavirus type 2; SPL: Sound pressure level; SRT: Speech reception threshold; STI: Speech transmission index

Supplementary Information

The online version contains supplementary material available at https://doi.org/10.1186/s12995-021-00334-y.

Additional file 1: Supplemental Fig. S1. Dampening function (i.e. frequency response) of the surgical face mask (light grey) and the FFP3 respirator (black).

Acknowledgements

The authors would like to thank Annette Lehn (Institute of Biostatistics and Mathematical Modelling, Goethe University, Frankfurt am Main) for her support in statistical data analysis and all subjects for their participation in this study.

Authors' contributions

RW, TS and MD conceived and designed the study. RW and TW² conducted audiological analyses and collected data. RW, LG, NF, TW² and MD performed the data analysis and visualization. TW², TW⁵ and TS provided resources for the study. RW and MD wrote the original draft of the manuscript. RW, LG, TW², NF, BH, KZ, TW⁵, SW, VAJK, SC, TS and MD made substantial contributions to data interpretation and/or review and editing of the manuscript. All authors read and approved the final manuscript.

Funding

Open Access funding enabled and organized by Projekt DEAL.

Availability of data and materials

The datasets used and/or analysed during the current study are available from the corresponding author on reasonable request.

Declarations

Ethics approval and consent to participate

The Goethe University Hospital Frankfurt's ethics committee reviewed and approved this study (No. 20–909). Written informed consent was obtained from all subjects (health care professionals) who volunteered to participate in this study.

Consent for publication

Written informed consent for publication of study data was obtained from all participants.

Competing interests

The authors declare that they have no competing interests.

Author details

¹Department of Otorhinolaryngology, Head and Neck Surgery, University Hospital Frankfurt, Goethe University, Theodor-Stern-Kai 7, 60590 Frankfurt/M, Germany. ²Division of Audiological Acoustics, Department of Otorhinolaryngology, Head and Neck Surgery, University Hospital Frankfurt, Goethe University, Theodor-Stern-Kai 7, 60590 Frankfurt/M, Germany. ³Institute of Biostatistics and Mathematical Modelling, Goethe University, Theodor-Stern-Kai 7, 60590 Frankfurt/M, Germany. ⁴Department of Anesthesiology, Intensive Care Medicine & Pain Therapy, University Hospital Frankfurt, Goethe University, Theodor-Stern-Kai 7, 60590 Frankfurt/M, Germany. ⁵Department of Internal Medicine II, Infectious Diseases, University Hospital Frankfurt, Goethe University, Theodor-Stern-Kai 7, 60590 Frankfurt/M, Germany. ⁶Occupational Health Service, University Hospital Frankfurt Goethe University, Theodor-Stern-Kai 7, 60590 Frankfurt/M, Germany. ⁷Institute for Medical Microbiology and Infection Control, University Hospital Frankfurt, Goethe University, Paul-Ehrlich-Str. 40, 60596 Frankfurt/M, Germany. ⁸University Center of Competence for Infection Control of the State of Hesse, Goethe University, Paul-Ehrlich-Str. 40, 60596 Frankfurt/M, Germany. ⁹Institute of Medical Virology, University Hospital Frankfurt, Goethe University, Paul-Ehrlich-Str. 40, 60596 Frankfurt/M, Germany. ¹⁰German Centre for Infection Research, External partner site Frankfurt, Frankfurt/M, Germany. ¹¹Fraunhofer Institute for Molecular Biology and Applied Ecology (IME), Branch Translational Medicine and Pharmacology, Frankfurt, Germany.

Received: 19 April 2021 Accepted: 21 September 2021 Published online: 30 September 2021

References

- Zhu N, Zhang D, Wang W, Li X, Yang B, Song J, et al. A novel coronavirus from patients with pneumonia in China, 2019. N Engl J Med. 2020;382(8): 727–33. https://doi.org/10.1056/NEJMoa2001017.
- Qiu C, Deng Z, Xiao Q, Shu Y, Deng Y, Wang H, et al. Transmission and clinical characteristics of coronavirus disease 2019 in 104 outside-Wuhan patients, China. J Med Virol. 2020;92(10):2027–35. https://doi.org/10.1002/jmv.25975.

- Huang C, Wang Y, Li X, Ren L, Zhao J, Hu Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020;395(10223):497–506. https://doi.org/10.1016/S0140-6736(2 0)30183-5.
- World Health Organization. Coronavirus disease 2019 (COVID-19) Situation Report - 73. 2020. https://www.who.int/docs/default-source/coronaviruse/ situation-reports/20200402-sitrep-73-covid-19.pdf. Accessed 14 Apr 2021.
- Patel ZM, Fernandez-Miranda J, Hwang PH, Nayak JV, Dodd RL, Sajjadi H, et al. Letter: precautions for endoscopic Transnasal Skull Base surgery during the COVID-19 pandemic. Neurosurgery. 2020;87(1):E66–7. https://doi.org/1 0.1093/neuros/nyaa125.
- Zou L, Ruan F, Huang M, Liang L, Huang H, Hong Z, et al. SARS-CoV-2 viral load in upper respiratory specimens of infected patients. N Engl J Med. 2020;382(12):1177–9. https://doi.org/10.1056/NEJMc2001737.
- Workman AD, Welling DB, Carter BS, Curry WT, Holbrook EH, Gray ST, et al. Endonasal instrumentation and aerosolization risk in the era of COVID-19: simulation, literature review, and proposed mitigation strategies. Int Forum Allergy Rhinol. 2020;10(7):798–805. https://doi.org/10.1002/alr.22577.
- Tran K, Cimon K, Severn M, Pessoa-Silva CL, Conly J. Aerosol generating procedures and risk of transmission of acute respiratory infections to healthcare workers: a systematic review. PLoS One. 2012;7(4):e35797. https:// doi.org/10.1371/journal.pone.0035797.
- Loth AG, Guderian DB, Haake B, Zacharowski K, Stöver T, Leinung M. Aerosol exposure during surgical tracheotomy in SARS-CoV-2 positive patients. Shock. 2020;55(4):472–8. https://doi.org/10.1097/SHK.000000000001655.
- Boccalatte LA, Larrañaga JJ, Perez Raffo GM, Teijido CA, García Fornari G, Staneloni MI, et al. Brief guideline for the prevention of COVID-19 infection in head and neck and otolaryngology surgeons. Am J Otolaryngol. 2020; 41(3):102484. https://doi.org/10.1016/j.amjoto.2020.102484.
- 11. Roberts V. To PAPR or not to PAPR? Can J Respir Ther. 2014;50(3):87-90.
- Patel ZM, Fernandez-Miranda J, Hwang PH, Nayak JV, Dodd RL, Sajjadi H, et al. In reply: precautions for endoscopic Transnasal Skull Base surgery during the COVID-19 pandemic. Neurosurgery. 2020;87(2):E162–3. https:// doi.org/10.1093/neuros/nyaa156.
- Givi B, Schiff BA, Chinn SB, Clayburgh D, Iyer NG, Jalisi S, et al. Safety recommendations for evaluation and surgery of the head and neck during the COVID-19 pandemic. JAMA Otolaryngol Head Neck Surg. 2020;146(6): 579–84. https://doi.org/10.1001/jamaoto.2020.0780.
- Canadian Society of Otolaryngology- Head & Neck Surgery. Guidance for Health Care Workers Performing Aerosol Generating Medical Procedures during the COVID-19 Pandemic Endorsed by the CSO-HNS Executive Committee. 2020. https://www.entcanada.org/wp-content/uploads/ Protocol-for-COVID-and-AGMP-3-iw-mailer.pdf. Accessed 14 Apr 2021.
- American Academy Of Otolaryngology- Head and Neck Surgery. Academy Supports CMS, Offers Specific Nasal Policy. 2020. https://www.entnet.org/ content/academy-supports-cms-offers-specific-nasal-policy-1. Accessed 14 Apr 2021.
- Kempfle JS, Panda A, Hottin M, Vinik K, Kozin ED, Ito CJ, et al. Effect of powered air-purifying respirators on speech recognition among health care workers. Otolaryngol Head Neck Surg. 2021;164(1):87–90. https://doi.org/1 0.1177/0194599820945685.
- Palmiero AJ, Symons D, Morgan JW, Shaffer RE. Speech intelligibility assessment of protective facemasks and air-purifying respirators. J Occup Environ Hyg. 2016;13(12):960–8. https://doi.org/10.1080/15459624.2016.1200723.
- HAHLBROCK KH. Speech audiometry and new word-tests. [Uber Sprachaudiometrie und neue Wörterteste]. Arch Ohren Nasen Kehlkopfheilkd. 1953;162(5):394–431. https://doi.org/10.1007/BF02105664.
- Wagener K, Brand T, Kühnel V, Kollmeier B. Entwicklung und Evaluation eines Satztests f
 ür die deutsche Sprache I–III: Design, Optimierung und Evaluation des Oldenburger Satztests. Z Audiol. 1999;38:4–15 44–56, 86–95.
- Seedorff M, Oleson J, McMurray B. Maybe maximal: good enough mixed models optimize power while controlling type I error; 2019. https://doi. org/10.31234/osf.io/xmhfr.
- van Doremalen N, Bushmaker T, Morris DH, Holbrook MG, Gamble A, Williamson BN, et al. Aerosol and surface stability of SARS-CoV-2 as compared with SARS-CoV-1. N Engl J Med. 2020;382(16):1564–7. https://doi. org/10.1056/NEJMc2004973.
- Stadnytskyi V, Bax CE, Bax A, Anfinrud P. The airborne lifetime of small speech droplets and their potential importance in SARS-CoV-2 transmission. Proc Natl Acad Sci U S A. 2020;117(22):11875–7. https://doi.org/10.1073/pna s.2006874117.

- Mick P, Murphy R. Aerosol-generating otolaryngology procedures and the need for enhanced PPE during the COVID-19 pandemic: a literature review. J Otolaryngol Head Neck Surg. 2020;49(1):29. https://doi.org/10.1186/s40463-020-00424-7.
- Howard BE. High-risk aerosol-generating procedures in COVID-19: respiratory protective equipment considerations. Otolaryngol Head Neck Surg. 2020;163(1):98–103. https://doi.org/10.1177/0194599820927335.
- Ing EB, Xu QA, Salimi A, Torun N. Physician deaths from corona virus (COVID-19) disease. Occup Med (Lond). 2020;70(5):370–4. https://doi.org/1 0.1093/occmed/kqaa088.
- Weiss R, Guchlerner L, Loth AG, Leinung M, Wicker S, Kempf VAJ, et al. Typical symptoms of common otorhinolaryngological diseases may mask a SARS-CoV-2 infection. Eur Arch Otorhinolaryngol. 2021;278(9):1–8. https:// doi.org/10.1007/s00405-021-06726-4.
- Chen JX, Workman AD, Chari DA, Jung DH, Kozin E, Lee DJ, et al. Demonstration and mitigation of aerosol and particle dispersion during mastoidectomy relevant to the COVID-19 era. Otol Neurotol. 2020;41(9): 1230–9. https://doi.org/10.1097/MAO.000000000002765.
- Chee VWT, Khoo ML-C, Lee SF, Lai YC, Chin NM. Infection control measures for operative procedures in severe acute respiratory syndrome-related patients. Anesthesiology. 2004;100(6):1394–8. https://doi.org/10.1097/ 00000542-200406000-00010.
- Wolf T, Kann G, Becker S, Stephan C, Brodt H-R, de Leuw P, et al. Severe Ebola virus disease with vascular leakage and multiorgan failure: treatment of a patient in intensive care. Lancet. 2015;385(9976):1428–35. https://doi. org/10.1016/S0140-6736(14)62384-9.
- European Centre for Disease Prevention and Control. Risk assessment: SARS-CoV-2 - increased circulation of variants of concern and vaccine rollout in the EU/EEA, 14th update. 2020. https://www.ecdc.europa.eu/en/publica tions-data/covid-19-risk-assessment-variants-vaccine-fourteenth-updatefebruary-2021. Accessed 14 Apr 2021.
- Sandaradura I, Goeman E, Pontivivo G, Fine E, Gray H, Kerr S, et al. A close shave? Performance of P2/N95 respirators in healthcare workers with facial hair: results of the BEARDS (BEnchmarking adequate respiratory DefenceS) study. J Hosp Infect. 2020;104(4):529–33. https://doi.org/10.1016/j.jhin.2020. 01.006.
- Lee HP, Wang DY. Objective assessment of increase in breathing resistance of N95 respirators on human subjects. Ann Occup Hyg. 2011;55:917–21. https://doi.org/10.1093/annhyg/mer065.
- Brown VA, van Engen KJ, Peelle JE. Face mask type affects audiovisual speech intelligibility and subjective listening effort in young and older adults. Cogn Res Princ Implic. 2021;6(1):49. https://doi.org/10.1186/s41235-021-00314-0.
- Toscano JC, Toscano CM. Effects of face masks on speech recognition in multi-talker babble noise. PLoS One. 2021;16(2):e0246842. https://doi.org/1 0.1371/journal.pone.0246842.

Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Ready to submit your research? Choose BMC and benefit from:

- fast, convenient online submission
- thorough peer review by experienced researchers in your field
- rapid publication on acceptance
- support for research data, including large and complex data types
- gold Open Access which fosters wider collaboration and increased citations
- maximum visibility for your research: over 100M website views per year

At BMC, research is always in progress.

Learn more biomedcentral.com/submissions

